



## Fractional Co2 Laser Beauty Machine 40W For Beauty Salon And Home

Our Product Introduction

For more products please visit us on [hifubeauty-machine.com](http://hifubeauty-machine.com)


### Basic Information

- Place of Origin: China Guangzhou
- Brand Name: LZR
- Certification: CE
- Model Number: T9
- Minimum Order Quantity: 1
- Price: confer
- Packaging Details: Standard export packing
- Delivery Time: 5-8days
- Payment Terms: L/C, D/A, D/P, T/T, Western Union, MoneyGram
- Supply Ability: 100 Piece/Pieces per Month



### Product Specification

- Name: Fractional Co2 Laser Beauty Machine
- Wavelength Range: 430/480/530/560/640/690 ~ 1200nm
- Usage: Beauty Salon, Spa, Home
- Pulse Width: 0.02~0.99s
- Screen: 8-inch Touchscreen
- Treatment Probes: 5 Pieces
- Energy Density: 1-10J/cm<sup>2</sup>
- Style: Portable
- Output Power: 40W
- Highlight: Fractional Co2 Laser Beauty Machine, Fractional Co2 Laser Machine, Fractional Co2 Laser Device 40W



### More Images





## Product Description

### 40W 8-Inch Touchscreen Fractional Co2 Beauty Salon Machine For Beauty Salon And Home

#### Product Description:

One of the most impressive features of this product is its wavelength range of 430/480/530/560/640/690 ~ 1200nm. This means that it can effectively target different skin types and hair colors, making it suitable for a wider range of users.

The Diode Laser Beauty Machine comes with 5 treatment probes and 5 treatment heads, each designed for specific areas of the body. This makes it easier to target specific areas and achieve better results.

Using this IPL Laser Beauty Machine is simple and convenient. It features a user-friendly interface that allows for easy navigation and control. The device is also easy to use, with simple instructions that can be followed by anyone.

With regular use, the Diode Laser Beauty Machine can provide long-lasting results. Users can expect to see a reduction in hair growth over time, resulting in smoother, silkier skin.

Overall, the Diode Laser Beauty Machine is a great investment for anyone looking for a home laser hair removal machine. Its portability, adjustable spot diameter, and range of treatment probes and heads make it a versatile and effective option for a wide range of users.

#### Features:

Product Name: Laser Beauty Machine

Model: LB-100

Output Power: 40W

Screen: 8-inch Touchscreen

Energy Density: 1-10J/cm<sup>2</sup>

Screen Size: 15.6 Inches

This Fractional Co2 Laser Beauty Machine is perfect for use in a professional spa or as a Home Beauty Machine. With an Output Power of 40W and an Energy Density of 1-10J/cm<sup>2</sup>, you can be sure that this machine will give you the results you desire. The 8-inch Touchscreen makes it easy to control and customize your treatments, while the large 15.6-inch screen size ensures that you can see everything clearly.

#### Technical Parameters:

Frequency	1-6Hz
Energy Density	1-10J/cm <sup>2</sup>
Laser Bar Number	3*755nm; 4*808nm; 3*1064nm
Screen	8-inch Touchscreen
Pulse Width	0.02~0.99s
Application	For Commercial
Spot Diameter	6mm-14mm Adjustable
Wavelength Range	430/480/530/560/640/690 ~ 1200nm
Treatment Probes	5 Pieces
Style	Portable

#### Applications:

Thanks to its portable design, the Hair Removal Laser Machine can be easily transported to various locations, making it an ideal option for mobile beauty professionals. With a high frequency of 1-6Hz, this machine can deliver rapid and efficient results for a wide range of applications.

One of the key benefits of the Hair Removal Laser Machine is its adjustable pulse width. With a range of 0.02~0.99s, you can customize your treatments to suit your specific needs and achieve the best possible results. Additionally, the machine offers a wavelength range of 430/480/530/560/640/690 ~ 1200nm, which can be used to target a variety of skin concerns.

When it comes to product application occasions and scenarios, the Hair Removal Laser Machine is suitable for a variety of settings. For professional use, it can be used in beauty salons and spas to provide a range of laser treatments to clients. These may include hair removal, skin rejuvenation, acne treatment, and more.

For those who prefer to enjoy the benefits of laser treatments at home, the Hair Removal Laser Machine is an excellent choice. It can be used to achieve smooth, hair-free skin in the comfort of your own home, without the need for expensive salon treatments. Additionally, the machine can be used to address a range of skin concerns, such as fine lines, wrinkles, and uneven skin tone.

Overall, the Hair Removal Laser Machine is a versatile and effective tool that can be used in a variety of settings. Whether you're a beauty professional or a beauty enthusiast, this device is an excellent investment that can help you achieve the results you desire.















Customize your Diode Laser Beauty Machine with our product customization services. Choose from the following options:

Screen: 8-inch Touchscreen

Treatment Probes: 5 Pieces

Output Power: 40W

Style: Portable

Spot Diameter: 6mm-14mm Adjustable

Whether you're looking for a Home Use Beauty Machine or an E Light Beauty Machine, our customization services can help you get the most out of your Laser Beauty Machine.



**Guangzhou Lizhirun Technology Development Co., Ltd**



18688905627



cenglinlin25@gmail.com



hifubauty-machine.com

Room 502, 5th Floor, Building F, Huixin Industrial Zone, No. 89, Hetai Road, Hebian Village, Baiyun District,  
Guangzhou